

Merilyzer | AutoQuant 100i™ / AutoQuant 100 Amara™

Fully Automated Biochemistry Analyzer

AutoQuant 100i™ & AutoQuant 100 Amara™ is a compact, automated random access analyzer with throughput of 100 tests / hour, ideal for small to mid sized laboratory.

Product code: FAAQ1-02

- Throughput: 100 Tests/hour
- Pre loaded test protocols
- Measurement Principles: Photometry (Rate & End point), Turbidimetry Immunoassay
- Random Access with STAT function
- 22 dedicated sample positions & 22 onboard cooled reagent positions
- **60 dismountable, easily replaceable and reusable cuvettes with automatic blank testing**
- 7-step automatic washing and automatic cuvette drying function
- **Halogen tungsten lamp with 9 channels optic fibres to enhance optical accuracy & life of interference filters**
- 9 wavelengths: 340, 405, 450, 510, 546, 578, 600, 670, 700 nm
- Calibration Points: K-Factor, Linear (one, two and multi point), Logit-log, Spline, Exponential, Polynomial
- **User defined Reflex Testing**
- **Low water consumption of 3.5 - 4.5 litres / hour**
- Operating System: Windows 10 IOT (2019)



SCAN FOR BROCHURE

