

# Merilyzer | AutoQuant 200i™ / AutoQuant 200 Excelus®

## Fully Automated Biochemistry Analyzer

AutoQuant 200i™ & AutoQuant 200 Excelus® is a fully automated, random access, patient prioritized analyzer with throughput of 200 tests / hour, making it ideal for small to mid sized laboratory.

Product code: FAAQ1-02

- Throughput: 200 Tests/hour
- Random Access with STAT function
- Measurement Principles: Photometry (Rate & End point), Turbidimetry Immunoassay
- 30 dedicated sample positions & 30 onboard cooled reagent positions
- **60 dismountable, easily replaceable and reusable cuvettes**
- **Halogen tungsten lamp with 9 channels optic fibres to enhance optical accuracy & life of interference filters**
- 9 wavelengths: 340, 405, 450,510, 546, 578, 600, 670, 700 nm
- Calibration Points: K-Factor, Linear (one, two and multi point), Logit-log, Spline, Exponential, Polynomial
- **User defined Reflex Testing**
- On board Laundry System, 7-step automatic washing
- **Low water consumption of 5-6 litres / hour**
- Operating System: Windows 10 IOT (2019)



SCAN FOR BROCHURE

