



ALL NEW RE-DESIGNED



HC-300 **NEW**
Urine Analyzer

DIRUI


DIRUI

Multi-scenario application

- o Beyond imagination

Portable

- o Easy to carry



Performance

- o A smart compact device

Comprehensive parameter: 10~14 chemistry parameters (Optional)
A/C (Calculation parameter)

Multi language compatibility: EN/RU/ES/IT/PT/DE/TR/FR/PL

User friendly design: Color touch screen



Technical Specification

Display	4.3" color touch screen
Throughput	60T/H or 120T/H
Test parameters	UBG, BIL, KET, CRE, BLD, PRO, MALB, NIT, LEU, GLU, SG, pH, VC, Ca
Principle	Photoelectric Colorimetry
Wavelengths	525nm, 572nm, 610nm and 660nm
Storage capacity	>5000, SD card result export
Languages	English, Spanish, Portuguese, Russian, Polish, Italian, Turkish, German, French
Dimensions	265×140×125mm
Power	12V, 110/220V power supply
Ambient Temperature	15°C~35°C
Relative Humidity	75% or less
Atmospheric Pressure	76kPa~106kPa

DIRUI

DIRUI INDUSTRIAL CO.,LTD.

3333 Yiju Street,New&High Tech. Development Zone Changchun,Jilin 130103,P.R.China
E-mail:marketing@dirui.com.cn Http://www.dirui.com.cn/en
Tel: +86(431)85083721

*Specifications subject to change without notice.
20210426

CE



Official WeChat
official account