





# FLAT PANEL DETECTOR

SEE BETTER, TREAT BETTER



## BEST OUTCOMES

- Optimal workflow for urology, orthopedics, and general surgery procedures
- Minimal dose per session
- Fast stone identification & diagnosis

# ENVISION HEALING

## HIGH IMAGE QUALITY

- Wide **dynamic range** with **enhanced contrast** and depth
- High visibility with low noise
- Low distortion and no aberrations

## Small Footprint

- 26.1cm x 25.6cm x 5.0cm
- Fits in tight spaces

## Large Field of View

- 9" x 9" (22.9 x 22.9cm)
- Square-shaped
- Clearly see anatomic details

## Low Dose Exposure

- Remote diagnosis
- Zero radiation to medical staff

## Dose Control

- Ability to adjust kV and mA; select manual or automatic dose control according to image feedback
- Foot switch controls dosage time

## Easy Positioning

- Autolocation feature of stone treatment area for exact planning & positioning
- Ergonomic handles for easy maneuvering

## Powerful

- 40 KHz continuous fluoroscopy frequency
- 3.5 KW nominal power



203 Perry Parkway, Suite #6,  
Gaithersburg, MD 20877, USA  
Tel: +1 (301) 944-1575  
Fax: +1 (301) 972-6098  
E-mail: [marketing@medispec.com](mailto:marketing@medispec.com)  
[www.medispec.com](http://www.medispec.com)

