


S150

AUTOMATIC
CHEMISTRY ANALYZER



ZHUHAI SENLONG BIOTECH CO.,LTD

Senlo[®]

Refrigerated Reagent Compartment



PRODUCT FEATURES

Mini Fully Auto Chemistry Analyzer
Discrete/random access
STAT sample priority

Easy to use
Internal computer with 13" displayer
All-in-one machine

Automatic Washing
No Washing (Option : SL150 B)

REACTION/REAGENT/SAMPLE HANDLING

40 cuvettes
25 sample positions+5 original tube
50 positions for reagent with cooling system

25 bottles with volume of 6ml
25 bottles with volume of 24ml
LED light source with 8 channel transmissometry

Liquid level detection / Collision protection

WORKING CONDITIONS

Temperature:10~30°C
Humidity: ≤85%

Power Supply: AC 100~240V 50/60HZ
Water Consumption: ≤5L/H

OTHER

Test speed : Constant 150T/H
Dimension : 352mm(L)*600mm(W)*385mm(H)
Net weight : 27.8KG

LCIWS

Online
Monitoring

Online
Maintenance



SENLO WEB

ZHUHAI SENLONG BIOTECH CO.,LTD

Tel: +86-756-3916978 Fax: +86-756-3882997 Customer Hotline: 400 6789 811 Web: www.senlo.com.cn

Address : Taichuan Industrial Park, 611 Huawei Road, Zhuhai

The technical information is for reference only. All the indexes are taking the technical requirements of the product as the criterion, and the right to modify the model and appearance of the product is reserved. Zhuhai Senlong biotechnology Co Ltd reserves right of final interpretation and decision. Version: SL2023